X 180 (NEXY 180) 175 CC

MOTORLANCE NEXY 180

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Contents of the resuscitation kit:

- 1- Medication bag
- 2- Syringes
- 3- Laryngoscope
- 4- 2-liter oxygen cylinder (with manometer, mask, and nasal cannula)
- 5- Sterile gauze
- 6- Ambu bag (with mouth opener)
- 7- Airway
- 8- Scissors
- 9- Various bandages (elastic and gauze/woven)
- 10- Various tapes (Leukoplast, non-woven/fixation tape, and hypoallergenic tape)
- 11- Latex gloves
- 12- Minor surgery set (forceps, tissue forceps, sharp-tipped scissors, scalpel blade, scalpel handle, scissors, Magill forceps, and kidney dish/receiver)
- 13- Manual suction device
- 14- Blood glucose meter
- 15- Fingertip pulse oximeter

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Safety and Maintenance Guidelines:

Pre-deployment Checks for the Motor-Ambulance:

Daily inspection and full preparation of the motor-ambulance include checking the equipment. This review ensures that everything is in its proper place and readily accessible.

To this end, you must check the following items on the motor-ambulance:

Fuel level, engine oil level, transmission fluid level, brake fluid, drive belt, horn, tire pressure, siren and alarms, and handheld radios. Additionally, perform a daily check to ensure all equipment is securely fastened and verify the status of medical supplies—including oxygen cylinders, resuscitation kits, dressings and bandages, cervical collars, other equipment, and emergency childbirth kits.

All battery-operated equipment—such as defibrillators, pulse oximeters, and laryngoscopes—must be tested daily. During this check, you determine whether the equipment is functioning correctly, whether the supplies are clean and sterile, and whether there is a sufficient quantity available.

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Preparation Phase:

Stock the motorcycle ambulance with medical supplies and equipment based on anticipated needs and requirements.

Prioritize items essential for patient care in life-threatening situations; these include equipment for airway management, artificial respiration, and oxygen delivery.

Include equipment for cardiac care, wound dressing, suturing, and blood pressure monitoring.

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Warranty and After-Sales Service:

Dear Customer: 

The registered trademark belongs exclusively to Nima P.A.S. Company. Products bearing this mark are covered by Nima P.A.S. after-sales services for the duration of the warranty period, calculated from the date of the sales invoice. Thanks to the high quality of the products and services provided, this trademark currently enjoys significant prestige in Iran.

Nima P.A.S. Company is equipped to address all customer-related inquiries in accordance with international standards, specifically the ISO 1002:2004 standard for customer complaint handling systems.

Steps to access Nima P.A.S. Company's after-sales services:

1- Contact Nima P.A.S. Company Customer Service:

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2- Scheduling the Service Appointment:

All requests are reviewed by specialists at the earliest opportunity. The customer will be contacted to schedule an appointment or arrange for the product to be sent to the head office or factory.

3- Actions Taken by the Technical Customer Service Unit:

Experienced specialists in the Technical Customer Service Unit will conduct a thorough inspection and troubleshooting of the product. They will promptly resolve the issue reported by the customer as well as any other potential defects. Subsequently, the product will undergo performance evaluations and specialized testing by the company's Quality Control Unit to ensure it is returned to the customer with its original quality standards.

4- Product Return:

The Customer Affairs Department will make the necessary arrangements to return the product to the customer as soon as possible.

5- Quality Assessment of Services Provided:

Your satisfaction with the quality and speed of the services provided by the After-Sales Service Department is a priority for Nima P.A.S. Company.

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Business Hours:

Dear Customers,

Nima P.A.S. Company is available to address all your inquiries from Saturday to Thursday, between 8:30 AM and 5:00 PM.

We welcome your comments, suggestions, and feedback.

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Sincerely, Nima P.A.S. Company (Trading, Equipping, and Custom Outfitting)

PASMED